Medical Devices

Rhino Clear Space

IT > MANUALE ISTRUZIONI D'USO
EN > OPERATING INSTRUCTIONS MANUAL
DE > BEDIENUNGSANLEITUNG
FR > MODE D'EMPLOI
NL > GEBRUIKSAANWIJZING
EL > ΕΓΧΕΙΡΙΔΙΟ ΟΔΗΓΙΩΝ ΧΡΗΣΗΣ
PL > INSTRUKCJA OBSŁUGI I KONSERWACJI
ES > MANUAL DE INSTRUCCIONES DE USO
SV > BRUKSANVISNING
HR > PRIRUČNIK S UPUTAMA ZA UPORABU

RC3-1

◀ MADE IN ITALY ▶

IT

EN

DE

FR

NL


EL

PL


ES

SV

HR

Schema di collegamento / Connection diagram /
Anschlussdiagramm / Schéma de branchement /
Verbindungsschema / Διάγραμμα σύνδεσης / Sche-
mat połączeń / Esquema de conexión / Anslut-
ningsschema / Dijagram povezivanja

A e B1 non sono in dotazione al dispositivo medico, fanno parte della dotazione originale del prodotto SpaceNeb già acquistata / A and B1 are not supplied with the medical device, they are part of the original set supplied with the previously purchased SpaceNeb product / A und B1 gehören nicht zum Lieferumfang des Medizinprodukts; sie sind Bestandteil der Originalausstattung des bereits erworbenen SpaceNeb-Produkts / A et B1 ne sont pas fournis avec le dispositif médical, ils font partie de l'équipement d'origine du produit SpaceNeb déjà acheté / A en B1 worden niet meegeleverd met het medische apparaat; ze maken deel uit van de originele

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Nasal douche

These user instructions are provided for device model RC3-1.

INTENDED PURPOSE. Medical Device for washing and treatment of the nasal cavities using saline solutions and/or medicines, the treatment and medicines must be prescribed by a doctor.

OPERATING INSTRUCTIONS. Treatment of upper airway disorders via the endonasal route (rhinitis, rhinosinusitis, sinusitis and/or as indicated by the treating physician). Physiological washing of the nasal cavities (with saline, isotonic, hypertonic solutions, or hot spring water). In children, it helps with nasal hygiene; in adults, it promotes physiological nasal rinsing and can be used before upper airway therapy to assist drug administration.

⚠️ CONTRAINDICATIONS. • The medical device should NOT be used for patients who are unable to breathe on their own or who are unconscious • Do not use the device in anaesthetic or assisted ventilation circuits.

INTENDED USERS. The device is intended to be used by legally qualified medical personnel/healthcare professionals (physicians, nurses, therapists, etc.). The device can be used directly by the patient.

⚠️ TARGET PATIENT GROUP. Adults, children of all ages, infants. Before using the device, the user manual must be read carefully and an adult responsible for safety must be present if the device is used by infants, children of any age or persons with limited abilities (e.g. physical, mental or sensory). It is the responsibility of the medical personnel to assess the patient's condition and capabilities in order to determine, when prescribing therapy, whether the patient is able to operate the device safely on their own or whether the therapy should be administered by a responsible person. Medical personnel should be responsible for assessing the use of the device on specific types of patients such as pregnant women, lactating women, infants, incapacitated persons or persons with limited physical capabilities.

OPERATING ENVIRONMENT. This device can be used in healthcare facilities, such as hospitals, outpatient clinics, etc., or even at home.

⚠️ WARNINGS ABOUT POSSIBLE MALFUNCTIONS. • Should your device fail to perform, please contact the authorised service centre for clarifications. • The manufacturer should be contacted to report problems and/or unexpected events relating to operation and if necessary for clarification on use and/or maintenance/cleaning. • Please also refer

EN

to the case history of faults and their resolution.

WARNINGS. • Use the device for therapeutic purposes only. This medical device is not intended as a life-saving device. Any other use is considered misuse and there may be hazards. The manufacturer is not liable for any improper use. • Always consult your general practitioner for identification of treatment. • Follow your doctor's or therapist's instructions regarding the possible use of types of medicine, dosage and indications for treatment. • Pay attention to the indications provided with the drug and avoid using the devices with substances and dilutions other than those recommended.

• If you should experience allergic reactions or other problems while using the device, stop using it immediately and consult your doctor. • Keep this manual safe for further reference. • If the packaging is damaged or opened, contact the distributor or service centre. • The time required to switch from storage to operating conditions is about 2 hours. • Repairs may only be carried out by personnel authorised by the manufacturer. Unauthorised repairs invalidate the warranty and may pose a hazard to the user. • The nasal douche must not be used as a toy. The utmost attention should be paid when used by children.

Risk of suffocation: Some components of the device are small enough to be swallowed by children; therefore, keep the device out of the reach of children.

Risk of ineffectiveness of the treatment: - In the case of substances that are too dense, dilution with a suitable saline solution may be necessary, as prescribed by a doctor. - Comply with the cleaning instructions for the nasal douche as described in the Hygienic preparation section. - Use only original Flaem spare parts; the manufacturer shall not be held liable if non-original spare parts are used.

Risk of infection: - Personal use of the nasal douche is recommended to avoid any risk of infection. - Follow hygienic preparation before each use. Ensure that the douche and its parts are not stored in the vicinity of other accessories or devices for different therapies (e.g. infusion devices). - Once therapy is complete, do not leave the used liquid (medication or saline solution) inside the spray chamber, and proceed with the hygienic preparation operations.

TROUBLESHOOTING

Problems	Solutions	
If the nasal douche delivers insufficiently or for a short time.	Check that the atomiser (C.3) has been positioned correctly. Make sure the nasal douche is clean; if not, clean the device again, carefully following the instructions for hygienic preparation.	
If the nasal douche does not spray.	Check that the solution has been correctly added into the spray chamber (C.4.1) of the ampoule body (C.4). Check that the amount of liquid does not exceed 10 ml. Check that the nasal douche is clean and, if it is not, perform cleaning again by carefully following the instructions for hygienic preparation.	
	The accessories are not properly connected to the compressor unit.	Check proper connection of the accessories to the compressor unit.
	The compressor unit delivers insufficient or no airflow.	Check that the compressor unit delivers airflow.

If, after checking the conditions described above, the device is still not working properly, we recommend contacting your trusted dealer or an authorised FLAEM service centre nearest to you. You can find a list of all Service Centres at <http://www.flaemnuova.it/it/info/assistenza>.

DISPOSAL. The nasal douche is to be disposed of as general waste after a sanitisation cycle.

Packaging

Product box

Heat-shrink film

REPORTING OF SERIOUS INCIDENTS. Please report any serious incidents occurring in connection with this device to the manufacturer and to the competent authority of the Member State where you reside. An event is considered serious if it causes or may cause, directly or indirectly, death or an unforeseen serious deterioration in a person's state of health.

INFORMATION ON RESTRICTIONS OR INCOMPATIBILITIES WITH CERTAIN SUBSTANCES.

• Interactions: The materials used in the device are biocompatible materials and comply with statutory regulations, however possible allergic reactions cannot be completely ruled out. • Use the medication as soon as possible once it has been opened and avoid leaving it in the nasal douche; once therapy has ended, do not leave the medication inside the nasal douche and proceed with hygienic preparation steps.

SYMBOLS ON DEVICE OR PACKAGING

 Unique device identifier	 Manufacturer	 Model number
 Type BF applied part	 Production date	 Batch code
 Attention	 Phthalate- and bisphenol-free	 See instructions for use
 Medical CE Marking ref. Regulation 2017/745 EU as amended	 Atmospheric pressure limits	 Moisture limits
 Medical device.	 Temperature limits.	

TECHNICAL FEATURES. Model: RC3-1

- Maximum spray chamber capacity: 10 ml.

- The nasal douche is used in combination with the Flaem appliance model P0905EM REF 1.

ENVIRONMENTAL CONDITIONS

Operating conditions:

Ambient temperature: Between +10°C and +40°C

Relative air humidity: Between 10% and 95%

Atmospheric pressure: Between 69 KPa and 106 KPa

Storage and transport conditions:

Ambient temperature: Between -25°C and +70°C

Relative air humidity: Between 10% and 95%

Atmospheric pressure: Between 69 KPa and 106 KPa

These indications refer to Aerosol therapy device model P0905EM REF 1. A4.1 - Speed level indicator LED	speed 1 	speed 2 	speed 3 
Dispensing time of 10 ml *	50"	1'	1' 10"
MMD [μm]**	>10		
* (Tests performed with saline solution 0.9% NaCl according to Flaem internal procedure) ** The MMD value shown was taken with the Malvern Spraytec instrument and refers to the use of saline (0.9% NaCl), so it may vary depending on the medication used. The value indicated does not apply to medicines dispensed in suspension or with high viscosity. In that case, information must be obtained from the supplier of the medication.			

- Dimensions (1): 5x5x11h cm.

- Weight (1): 40 g.

- Type BF applied parts are the 3 nasal adaptors (C.1A, C.1B, C.1C).

DURATION. Model: RC3-1. The service life of the device is one year

EQUIPMENT AND INFORMATION ON MATERIALS

C.1A - Nasal adaptor 'A', suitable for children up to 3 years of age C.1B - Nasal adaptor 'B', suitable for 3 to 12 years of age C.1C - Nasal adaptor 'C', suitable for over 12 years of age	Polypropylene + thermo-plastic rubber
C.2 - Separator (Recovery chamber) C.3 - Atomiser C.4 - Douche body with ergonomic handle C.4.1 - 10 ml spray chamber C.4.2 - Air inlet	Polypropylene
IMPORTANT NOTE: Remove the identification label on the packaging and affix it in the appropriate spaces on page 2.	

INSTRUCTIONS FOR USE

Before each use, wash your hands thoroughly and clean the NASAL DOUCHE as described in the Hygienic preparation section. During application, it is recommended to adequately protect yourself from any dripping.

HOW TO USE IN COMBINATION WITH AEROSOL THERAPY DEVICE MODEL P0905EM REF 1

Mode 1 - Direct use (Appliance held by the user) (P. 4). The device is held directly by the user during the delivery of the medication or saline solution.

Mode 2 - Use with connecting hose B1 (P. 4). The compressor is placed on a flat, stable surface (e.g. a table). The nasal douche is connected to the main unit using the dedicated flexible tube, allowing the nasal douche to be brought close to the patient's nose. This configuration is especially ideal for: paediatric patients, elderly or dependent patients who may have difficulty holding the device, or situations in which greater freedom of movement is required.

1. Insert the atomiser (C.3) into the spray chamber (C.4.1).
2. Pour the solution (max 10 ml) into the spray chamber (C.4.1).
3. Insert the separator (C.2) into the douche body (C.4) and rotate it clockwise to secure it.
4. Then insert the age-specific adaptor or the most suitable one for the size of your nostrils (C.1A/C.1B/C.1C).
5. **Mode 1.** Insert the assembled douche (C) into the compressor unit; the symbol (P) must be positioned as in the image (P. 4).

Mode 2. Connect the assembled douche (C) to the compressor unit using the hose (P. 3).
6. Select the desired speed by pressing button A4 depending on the type of treatment you wish to carry out, see section **"TECHNICAL FEATURES. Model: RC3-1"**.

7. Bring the nasal douche close to your nose by placing the adaptor hole against one nostril and turn on the appliance by pressing the switch (A1).

8. Inhale and exhale.

9. During treatment, it is preferable to keep the other untreated nostril closed with a finger.

10. Repeat from step 6. to step 9. for the other nostril as well, distributing the solution evenly between the two nostrils.

11. In the case of children or impaired persons, and however unable to co-ordinate the various steps in 10., it is recommended to carry out the application alternating between the two nostrils every 20 seconds approximately, with the necessary pauses to expel the

mucus.

12. What is described in the previous points is useful for nasal douching using saline solutions. It is advisable to carry out the drug treatment after initial washing with saline solution to achieve a more effective distribution of the drug on the nasal mucous membranes.

13. When application is over, switch off the appliance and carry out cleaning following the instructions provided in section “**HYGIENIC PREPARATION**”.

HYGIENIC PREPARATION

Disassemble the Nasal Douche as shown in the connection diagram.

Then proceed according to the instructions below.

Sanitisation. Before and after each use, sanitise the components, choosing one of the methods described below.

method A: sanitise the components under warm drinking water (approx. 40°C) with mild dish-washing detergent (non-abrasive).

method B: sanitise the components in a dishwasher on a hot cycle.

method C: sanitise the components by immersing them in a 50% water and 50% white vinegar solution, then rinse thoroughly with warm drinking water (approx. 40°C).

If you also wish to carry out cleaning by DISINFECTION, skip to the DISINFECTION section.

After sanitising the components, shake them vigorously and place them on a paper towel, or alternatively dry them with a hot air jet (e.g. hair dryer).

Disinfection. After sanitising the components, disinfect them using one of the methods described below.

method A: obtain a disinfectant of the electrolytic chloride disinfectant (active ingredient: sodium hypochlorite), specifically for disinfection, available from pharmacies.

Execution: - Fill a container big enough to hold all of the parts being disinfected with a solution of potable water and disinfectant, according to the proportions indicated on the packaging of the disinfectant. - Completely immerse each part in the solution, taking care to avoid the formation of air bubbles on the parts. Leave the components immersed for the period of time indicated on the disinfectant packaging, and associated with the concentration chosen for the preparation of the solution. - Remove the disinfected parts and rinse abundantly with warm potable water. - Dispose of the solution according to the instructions provided by the disinfectant manufacturer.

method B: disinfect the components by boiling them in water for 10 minutes; use demineralised or distilled water to avoid lime deposits.

method C: disinfect the components using a steam baby bottle steriliser (not a microwave type). Carry out the process strictly according to the hot bottle steriliser's instructions. For disinfection to be effective, choose a steriliser with an operating cycle of at least 6 minutes. After disinfecting the components, shake them vigorously and place them on a paper towel, or alternatively dry them with a hot air jet (e.g. hair dryer).

Sterilisation

Sterilisation for clinical or hospital use is the responsibility of the facility performing the sterilisation and must be carried out according to the instructions in this manual.

The sterilisation procedure described in this section is only effective on treated components if all the points are adhered to and only if the components to be treated are sanitised beforehand, and is validated in accordance with ISO 17665.

Equipment: Steam steriliser with fractionated vacuum and overpressure compliant with EN 13060.

Execution: Package each individual component to be treated in sterile barrier system or packaging in accordance with EN 11607-1. Place the packaged components in the steam steriliser. Perform the sterilisation cycle in accordance with the instructions for use of the equipment by selecting a temperature of 134°C and a time of 4 minutes first.

Storage: Store sterilised components according to the instructions for use of the system or sterile barrier packaging chosen.

**ADDITIONAL INFORMATION FOR CLEANING the spray chamber**

In order to prevent deposits from clogging the atomiser and spray chamber and thereby affecting the function of the nasal douche, it is recommended that cleaning and sanitising be carried out after each use, following the directions in the instructions for use supplied with the nasal douche. After thoroughly rinsing the parts under lukewarm tap water, shake them to remove any remaining drops of water. Then place the parts on a paper towel or dry them with a jet of hot air. **ATTENTION:** If your nasal douche does not atomise correctly after cleaning and sanitising, do not use sharp tools (e.g. screwdrivers, scissors, needles, hairpins, etc.) to try to clean or unclog the holes in the atomiser or spray chamber, as this may damage these parts and impair the proper functioning of the unit. Instead, perform additional cleaning and sanitisation cycles as described in the instruction manual and recommended above, until the nasal douche starts working properly again.

Table of intended methods / components

Components	C.1A	C.1B	C.1C	C.2	C.3	C4
HYGIENIC PREPARATION IN THE HOME						
Sanitisation						
method A	✓	✓	✓	✓	✓	✓
method B	✓	✓	✓	✓	✓	✓
method C	✓	✓	✓	✓	✓	✓
Disinfection						
method A	300	300	300	300	300	300
method B	40	40	40	40	40	40
method C	40	40	40	40	40	40
HYGIENIC PREPARATION IN A CLINIC OR HOSPITAL						
Sterilisation	40	40	40	40	40	40
✓: intended ✗: not intended	300: ✓MAX 300 TIMES			40: ✓MAX 40 TIMES		

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© 2026 FLAEM NUOVA®
FLAEM NUOVA S.p.A. non risponde di ulteriori estensioni del periodo di garanzia assicurate da terzi.
FLAEM NUOVA S.p.A.