
Neo Spire[®]

IT > MANUALE ISTRUZIONI D'USO
EN > USER MANUAL
DE > BEDIENUNGSANLEITUNG
FR > MODE D'EMPLOI
NL > GEBRUIKERSHANDLEIDING
PL > INSTRUKCJA OBSŁUGI
ES > MANUAL DE INSTRUCCIONES
EL > ΕΓΧΕΙΡΙΔΙΟ ΟΔΗΓΙΩΝ
TR > KULLANICI KILAVUZU

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Respiratory rehabilitation physiotherapy apparatus

These user instructions are provided for the device model P0713EM.

INTENDED PURPOSE. Medical device for the removal of tracheobronchial secretions and for respiratory physiotherapy, the treatment of which must be prescribed by a doctor.

OPERATING INSTRUCTIONS. Treatment of conditions of the respiratory system, in particular for the removal of tracheobronchial secretions (COPD, cystic fibrosis, bronchiectasis and/or as indicated by the treating physician).

⚠️ CONTRAINDICATIONS. • The medical device must NOT be used for patients who are unable to breathe on their own or who are unconscious. • Do not use the device in anaesthetic or assisted ventilation circuits. • The device must not be used for patients with pronounced arterial hypotension. The device must not be used on neonates.

INTENDED USERS. The devices are intended for use by legally authorised medical personnel/healthcare professionals (doctors, nurses, therapists, etc.). The device can be used directly by the patient for home therapy, but under the guidance of medical personnel or the legally authorised healthcare professional.

⚠️ TARGET PATIENT GROUP. Adults, children of all ages. Before using the device, the user manual must be read carefully and an adult responsible for safety must be present if the device is used by children of any age or persons with limited abilities (e.g., physical, mental or sensory). It is up to the medical personnel to assess the patient's condition and capabilities in order to determine, when prescribing therapy, whether the patient is able to operate the device safely on their own or whether the therapy should be administered by a responsible person.

Please refer to medical personnel for the evaluation of the use of the device on particular types of patients such as pregnant women, lactating women, incapacitated subjects or persons with limited physical capabilities.

OPERATING ENVIRONMENT. This device can be used in healthcare facilities, such as hospitals, outpatient clinics, etc., or even at home.

⚠️ WARNINGS REGARDING POSSIBLE MALFUNCTIONS. • Should your apparatus fail to perform, contact the authorised service centre for clarification. • The manufacturer should be contacted to report problems and/or unexpected occurrences relating to operation and if necessary for clarification of use and/or maintenance/hygienic preparation. • Please also refer to the case history of faults and their resolution.

WARNINGS. • Only use the device for its intended purpose. This medical device is not intended as a life-saving device. Any other use is considered misuse and there may be hazards. The manufacturer is not liable for any misuse. • Always consult your treating physician when determining the dosage of the treatment. • If you should experience allergic reactions or other problems while using the device, discontinue its use immediately and consult your treating physician. • Carefully store this manual in a safe place for further reference. • If the packaging is damaged or opened, contact your distributor or service centre. • Do not expose the apparatus to particularly extreme temperatures. • The time required to change from storage to operating conditions is approx. 2 hours. • It is forbidden to open and modify the device in any way. Repairs may only be carried out by personnel authorised by the manufacturer. Unauthorised repairs invalidate the warranty and may pose a hazard to the user.

• The apparatus must be transported inside the supplied bag and enclosed in its protective pouch.

• **Risk of suffocation:** Some components of the apparatus are small enough to be swallowed by children, so keep the apparatus out of the reach of children.

• **Risk of strangulation:** Do not use the supplied connecting hose and cables outside of their intended use, they could cause a strangulation hazard, take special care for children and people with special difficulties, often these people are not able to correctly assess the dangers.

• **Risk of fire:** This apparatus is not suitable for use in the presence of an anaesthetic mixture that is flammable with air, or with oxygen or nitrous oxide.

- Never obstruct the ventilation slots located on both sides of the device.

• **Risk of electrocution:** - Before first use, and periodically during the life of the product, check the integrity of the structure of the apparatus and the power cord to ensure that there is no damage; if it is damaged, do not plug it in and immediately take the product to an authorised service centre or your dealer.

- Keep the power cord away from animals (e.g. rodents), otherwise such animals could damage the insulation of the power cord. - Always keep the power cord away from hot surfaces. - Do not handle the plug and power cord with wet hands.

- Do not handle the apparatus with wet hands. Do not use the apparatus in humid environments (e.g.

while bathing or showering). Do not immerse the apparatus in water; if this happens, pull the plug immediately. Do not pull out or touch the apparatus while immersed in water, unplug it first. Take it immediately to an authorised FLAEM service centre or your dealer. - No maintenance/cleaning must be carried out during use.

Risk of ineffective treatment: - Use only original Flaem spare parts, no liability is accepted if non-original spare parts are used. - Observe hygienic preparation as described in the section Hygienic preparation. - Only use the apparatus in a dust-free environment, otherwise the therapy may be impaired. - Do not obstruct or insert objects into the filter and its housing in the apparatus.

Risk of infection: - Follow the hygienic preparation operations before and after each use. Ensure that the parts are not stored in the vicinity of other accessories or devices for different treatments (e.g. infusion devices). - The supplied parts are single-patient.

Risk of injury: - Do not place the apparatus on a soft supporting surface, such as a sofa, bed or tablecloth. - Always operate it on a hard, unobstructed surface.

WARNINGS ON RISKS OF INTERFERENCE DURING USE IN DIAGNOSTIC INVESTIGATIONS. This apparatus is designed to meet the current requirements for electromagnetic compatibility. As far as EMC requirements are concerned, electromedical devices require special care when being installed and used, and are therefore required to be installed and/or used in accordance with the manufacturer's specifications. Risk of potential electromagnetic interference with other devices. Mobile or portable RF radio and telecommunication devices (mobile phones or wireless connections) may interfere with the operation of electromedical devices. The apparatus may be susceptible to electromagnetic interference in the presence of other devices used for specific diagnosis or treatment. For more information visit www.flaemnuova.it.

TROUBLESHOOTING

Before performing any operation, switch off the apparatus and disconnect it from the mains.

Problem	Cause	Remedy
The apparatus does not switch on	Power cable not correctly plugged into the apparatus socket or mains power socket	Plug the power cable correctly into the sockets
The apparatus is not working	The dispenser has not been assembled correctly.	Disassemble and reassemble the dispenser correctly according to the connection diagram on the cover.
	The air hose is not properly connected to the apparatus	Check the correct connection between the air inlet of the apparatus and the accompanying accessories (see connection diagram on the cover).
	The air hose is bent or damaged or twisted	Unwind the hose and check it for crushing or punctures. Replace it if necessary.
	The air filter is dirty	Replace the filter
The apparatus is louder than usual	Filter cover not properly seated	Insert the Filter Cover properly into the housing

If, after checking the conditions described above, the apparatus still is not working properly, we recommend that you contact your trusted dealer or an authorised FLAEM service centre nearest to you. You can find a list of all Service Centres at <http://www.flaemnuova.it/it/info/assistenza>

DISPOSAL. Compressor unit

In accordance with Directive 2012/19/EC, the symbol on the appliance indicates that the apparatus to be disposed of (excluding accompanying components) is considered as waste and must therefore be subject to "separate collection". Therefore, the user must deliver this waste (or have it delivered) to the separate collection centres set up by the local authorities, or hand it over to the retailer when purchasing a new appliance of an equivalent type. Separate waste collection and subsequent treatment, recovery and disposal operations promote the production of equipment from recycled materials and limit the negative environmental and health effects caused by improper waste management. Unauthorised disposal of the product by the user entails the application of the administrative sanctions provided for in the laws transposing Directive 2012/19/EC of the member state or country where the product is disposed of.

Packaging

Product box

Tube packaging pouch

Heat-shrinkable nebuliser film and accompanying components

Sack packaging bag

REPORTING SERIOUS EVENTS. Serious events occurring in connection with this product must be reported immediately to the manufacturer or the competent authority. An event is considered serious if it causes or may cause, directly or indirectly, death or an unforeseen serious deterioration in a person's state of health.

SYMBOLS ON DEVICE OR PACKAGING			
Unique device identifier	Serial number of the apparatus	Class II device	Manufacturer
Type BF applied part	Production date	Batch code	Alternating current
Attention	Phthalate- and bisphenol-free	Keep dry	Model number
Medical device	Temperature limits	Moisture limits	Atmospheric pressure limits
Switched "ON"	When switching off the unit, the switch stops compressor operation on only one of the two power phases.		Before use: Caution check instructions for use
Switched "OFF"			Single patient
IP21	Protection rating of the enclosure: IP21. (Protected against solid bodies larger than 12 mm. Protected against access with a finger; Protected against vertically falling drops of water).		
CE 0051	Medical CE marking ref. regulation 2017/745 EU as amended		

INFORMATION ON RESTRICTIONS OR INCOMPATIBILITIES WITH CERTAIN SUBSTANCES.

• Interactions: The materials used in the device are biocompatible materials and comply with statutory regulations, however possible allergic reactions cannot be completely excluded.

TECHNICAL FEATURES. Model: P0713EM

Power supply: 230V ~ 50Hz 210VA

Air flow to compressor: 17 l/min approx

Noise level (at 1 m): 56 dB (A) approx

Operation: Continued

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Weight: 2.2 Kg

APPLIED PARTS

Type BF applied parts are:

accompanying component (B2.3)

ENVIRONMENTAL CONDITIONS

Operating conditions:

Ambient temperature: Between +5°C and +40°C

Relative air humidity: Between 10% and 95%

Atmospheric pressure: Between 69 KPa and 106 KPa

Storage and transport conditions:

Ambient temperature: Between -25°C and +70°C

Relative air humidity: Between 10% and 95%

Atmospheric pressure: Between 69 KPa and 106 KPa

DURATION

Main unit (A) Expected lifetime life 2000 hours.

For single-patient accompanying components (B, C1, C2) The average life expectancy is 1 year, however, it is advisable to replace them every 6 months in intensive use to ensure maximum therapeutic efficacy at all times.

APPLIANCE AND MATERIAL INFORMATION		
The equipment includes:		Information on materials
A -	Physiotherapy apparatus (main unit) A1 - Switch, A2 - Air inlet, A3 - Filter Cover A4 - Dispenser support, A5 - Carrying handle A6 - Socket for power cord A7 - Power cord	
	Accompanying components (single-patient) B1 -Connecting hose (main unit / dispenser) B2 - Dispenser B2.1 - Dispenser Body, B2.2 - Fitting, B2.3 - Mouthpiece	Polypropylene
D -	Comfortable and spacious carry bag	

OPERATING INSTRUCTIONS

Before each use, wash your hands thoroughly and clean your appliance as described in the section "HYGIENE PREPARATION". During application, it is recommended to adequately protect oneself from any dripping. The dispenser and accompanying components are for personal use only to avoid any risk of infection.

1. Plug the power cord (A7) into the socket (A6) and then connect it to a mains socket corresponding to the voltage of the apparatus. It must be positioned in such a way that disconnection from the power grid is not difficult.
2. Connect one end of the connecting hose (B1) to the apparatus. Connect the other end of the hose (B1) to the dispenser (B2).
3. Connect the accompanying components as shown in the "Wiring diagram" on the cover. Sit comfortably while holding the dispenser in your hand and place the mouthpiece to your mouth (Fig. 1).
4. Switch on the device by operating the switch (A1) and inhale and exhale deeply.
5. When the application is finished, switch off the apparatus and unplug it. The model P0713EM can be used in sessions of about 15 minutes several times during the day as prescribed by the treating physician.
6. You can insert the dispenser in the dispenser support (A4) if a break is needed during the therapy session; you will thus avoid contaminating the mouthpiece by placing it on dirty surfaces (Fig. 2).

ATTENTION: If, after the therapy session, an obvious deposit of moisture forms inside the tube (B1), detach the tube from the nebuliser and dry it with the ventilation of the compressor itself; this action prevents possible mould growth inside the tube.

To facilitate the connection of the connection tube (B1) to the compressor unit, act on the end of it by simultaneously rotating and inserting it, and for removal by rotating and extracting it.

HYGIENIC PREPARATION

Switch off the appliance before each hygienic preparation operation and disconnect it from the mains.

Apparatus (A) and outside of hose (B1). Use only a cloth moistened with antibacterial detergent (non-abrasive and free of solvents of any kind).

Accompanying components

The use of the dispenser and accompanying components must be personal to avoid any risk of infection. Dismantle the dispenser as shown in "Connection diagram" sect. B.

Then proceed according to the instructions below.

Sanitisation. Before and after each use, sanitise the accompanying components by choosing one of the methods provided in the table and described below.

method A: Sanitise the accompanying components under warm drinking water (approx. 40°C) with mild washing up liquid (non-abrasive).

method B: Sanitise the accompanying components in the dishwasher with the hot cycle (70°C).

method C: Sanitise the accompanying components by soaking in a solution of 50% water and 50% white vinegar, then rinse thoroughly with warm drinking water (approx. 40°C).

After sanitising the accompanying components, shake them vigorously and place them on a paper towel, or alternatively dry them with a hot air jet (e.g. hair dryer).

Disinfection. After sanitising the accompanying components, disinfect them by choosing one of the methods provided in the table and described below. Each method can be performed for a limited number of times (see figure in table).

method A: Obtain a disinfectant of the electrolytic chloride type (active ingredient: sodium hypochlorite), specifically for disinfection, available in all pharmacies.

Execution: - Fill a container of a suitable size to hold all the individual components to be disinfected with a solution of drinking water and disinfectant, observing the proportions indicated on the disinfectant packaging. - Completely immerse each individual component in the solution, taking care to avoid the formation of air bubbles in contact with the components. Leave the components submerged for the period of time indicated on the disinfectant packaging and associated with the concentration chosen for the preparation of the solution. - Retrieve the disinfected components and rinse them thoroughly with lukewarm drinking water. - Throw out the solution in accordance with the instructions provided by the disinfectant's producer.

method B: Disinfect the accompanying components by boiling them in water for 10 minutes; use demineralised or distilled water to prevent hard water deposits.

method C: Disinfect the accompanying components in a hot bottle steamer (not microwave). Carry out the process strictly according to the hot bottle steamer's instructions. For disinfection to be effective, choose a hot bottle steamer with an operating cycle of at least 6 minutes.

If you also want to carry out sterilisation, skip to the **Sterilisation** section.

After disinfecting the accompanying components, shake them vigorously and place them on a paper towel, or alternatively dry them with a hot air jet (e.g. hair dryer).

At the end of each use, store the apparatus complete with its accompanying accessories in a dry and dust-free place.

Sterilisation. Appliance: Steam steriliser with fractionated vacuum and overpressure compliant with EN 13060.

Execution: Package each individual component to be treated in sterile barrier system or packaging in accordance with EN 11607. Place the packaged components in the steam steriliser. Perform the sterilisation cycle in accordance with the instructions for use of the equipment by selecting a temperature of 134°C and a time of 4 minutes first.

Storage: Store the sterilised components according to the instructions for use of the system or sterile barrier packaging chosen.

The sterilisation procedure was validated in accordance with ISO 17665-1.

At the end of each use, store the apparatus complete with its accompanying accessories in a dry and dust-free place.

Table of planned methods / patient accessories			
Patient accessory	B2.3	B2.2	B2.1
HYGIENIC PREPARATION AT HOME			
Sanitisation			
method A	✓	✓	✓
method B	✓	✓	✓
method C	✓	✓	✓
Disinfection			
method A	300	300	300
method B	300	300	300
method C	300	300	300
HYGIENIC PREPARATION IN THE CLINIC OR HOSPITAL			
Disinfection			
method A	300	300	300
method B	300	300	300
Sterilisation	300	300	300
✓: envisaged \: not envisaged 300 :MAX 300 TIMES			

AIR FILTRATION

The apparatus is equipped with a suction filter, located under the Filter Cover (A3), which is to be replaced when it is dirty or changes colour. For replacement, lift the Filter Cover (A3) and pull out the Filter. Do not wash or re-use the same filter. Regular replacement of the filter is necessary to help ensure proper compressor performance. The filter must be checked on a regular basis. Contact your dealer or authorised service centre for replacement filters.

To replace the filter, pull it out as shown in (Fig.3).

The filter is designed so that it is always fitted in its housing. Do not replace the filter during use. **Only use original Flaem spare parts, no liability is accepted if non-original spare parts are used.**

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© 2025 FLAEM NUOVA®
FLAEM NUOVA S.p.A. non risponde di ulteriori estensioni del periodo di garanzia assicurate da terzi.
FLAEM NUOVA S.p.A.